

# NU SKIN® FACIAL SPA AND NU SKIN® CONDUCTIVE GEL

UNLOCK THE FUTURE OF SKIN CARE WITH  
NU SKIN'S LUXURIOUS IN-HOME SPA



## PRODUCT OVERVIEW

Nu Skin has unlocked the future of skin care with its in-home spa, the Nu Skin Facial Spa with Conductive Gel. Using microcurrent technology, this powerful combination stimulates and tones the skin resulting in an improved appearance.

Now you can easily perform three regimens per week and enjoy the benefits of a spa visit as a part of your regular anti-aging skin care routine—all in the comfort of your own home. With two interchangeable conductors, Facial Spa features settings for a 10-minute, full-face session and three five-minute, targeted area sessions to suit your needs, preferences, and lifestyle. Included Conductive Gel conducts the current and helps the Facial Spa device glide smoothly on the skin.

Nu Skin Facial Spa also includes a large, bright display and enhanced ergonomics, making it an easy-to-use, effective in-home spa system.

Indulge in a luxurious 10-minute spa experience for a more radiant, youthful appearance that will leave people wondering about your age—and your secret.

*Note: Please remember to read all directions before using the Facial Spa.*

## BENEFITS

The Facial Spa instrument is intended for skin stimulation for facial toning and improved appearance. The Facial Spa features microcurrent technology that works with the Nu Skin Conductive Gel.

- Stimulates skin for improved appearance
- Provides facial toning—improving facial contour and making the skin look more voluminous.
- Features microcurrent technology

## KEY FEATURES

- **Ergonomic design**—the conductors are specially designed to follow the contours of the face, while keeping the Conductive Gel in constant contact with your skin.
- **Lighted display**—an easy-to-read LCD display to inform you while the instrument is in active mode.
- **Audio indications**—the device beeps every 10 seconds and has “end treatment” and “low battery” notifications.
- **Large Conductor**—the large, contoured, and grooved surface area is designed for use during the full-face regimen or the forehead regimen.
- **Small Conductor**—the flat, smooth, oval shape centralizes surface contact area, helping deliver current to specific areas during the targeted area regimens.

## FACIAL SPA COMPATIBLE NU SKIN CONDUCTIVE GEL

Specifically formulated for use with the Nu Skin Facial Spa, Nu Skin Conductive Gel ensures adequate conductivity throughout your spa session, while helping the Nu Skin Facial Spa glide across the skin.

**Ingredients:** Water (Aqua), Glycerin, Propylene Glycol, Carbomer, Sodium Hydroxide, Phenoxyethanol, Methylparaben, Butylparaben, Ethylparaben, Isobutylparaben, Propylparaben

## FREQUENTLY ASKED QUESTIONS

**How long will the batteries last?**

In-house testing showed that batteries in the Facial Spa instrument lasted several months with normal usage.

# NU SKIN® FACIAL SPA AND NU SKIN® CONDUCTIVE GEL

## Why does the Facial Spa instrument beep?

When you begin a session, the device will beep indicating that your spa session is beginning. You will also hear beeps at 10-second intervals during your session indicating that the unit is on and working. The audible beeps are the **ONLY** indicator that the device is working. Once the spa session time reaches 30 seconds remaining, the spa will beep every five seconds until the spa session is complete, at which point the “spa session complete” alert will sound.

## How do I change the conductor attachments?

**To release a conductor:** With your free hand on the conductor being released, hold the instrument with the chrome activation panel facing you and press on the release button firmly with your thumb until the conductor is released from the instrument.

**To attach a conductor:** Hold the instrument with the display side up and hold the conductor with the positioning indicator (small etch) facing you. Press firmly to snap the conductor into place. If the conductor fails to attach, make certain the conductor is positioned correctly. Forcing the conductor could damage both the instrument and the conductor. The instrument will not work unless the conductor is attached correctly.

## Is Nu Skin Facial Spa a medical device?

Nu Skin has received FDA clearance to market the Nu Skin Facial Spa as a class II medical device for over-the-counter cosmetic use. The Facial Spa with its accompanying Conductive Gel is a facial toning device designed to stimulate and tone the skin for an improved appearance.

## How often should I use the Nu Skin Facial Spa?

You can conduct a Nu Skin Facial Spa regimen up to three times each week. Each regimen consists of two sessions—Session A and Session B. Facial Spa should be used as directed and only used with Conductive Gel.

## What are the session times?

Nu Skin Facial Spa features a 10-minute, full-face, basic regimen and three five-minute area regimens. The basic regimen has a five-minute Spa Session A and a five-minute Spa Session B. All of the area regimens include a two-minute Spa Session A and a three-minute Spa Session B.

## Can Facial Spa be used on the neck?

Nu Skin Facial Spa is registered for facial use only and should not be used on other parts of the body including the neck.

## To reduce the risk to your health: use product only for its intended use as described in the Facial Spa User Manual.

- If you are in the care of a physician, consult with your physician before using this device.
- Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Do not apply stimulation across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.
- Do not apply stimulation over painful areas. If you have painful areas, you should consult with your physician before using this device.
- Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- Do not apply stimulation over, or in proximity to, cancerous lesions or diseased tissues.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Do not apply stimulation while sleeping.
- Do not change, release, or attach conductors when the device is turned ON. You risk exposing yourself to the exposed electrical lead, which may result in electrical shock or burning of the skin.
- Do not conduct more than three Facia Spa regimens per week.
- Do not apply stimulation while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury.
- Adult use only—do not use the device on children; it has not been evaluated for pediatric use.
- Consult with your physician before using this device. The device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Apply stimulation only to normal, intact, clean, healthy skin.
- Do not operate while bathing or showering. To reduce the risk of shock, burn, fire, or injury, do not immerse your device in water. Do not place or store where it may fall or be pulled into a bathtub, shower, sink, or toilet.
- Do not apply stimulation on skin without Nu Skin® Conductive Gel. Using the device without Conductive Gel may burn your skin. Only use with Nu Skin® Conductive Gel.
- Have a mental implant or electronic implanted devices.
- Have a suspected or diagnosed epilepsy or suffer from seizures.
- Have metal braces on your teeth
- Are allergic to gels or metal or are sensitive to contact with chrome.
- Are pregnant or attempting to become pregnant.
- Have open wounds or sores.
- Are prone to rosacea.
- Have problematic or sensitive skin.

The Facial Spa and Nu Skin Conductive Gel are contraindicated for use in the following areas:

- Eyes, eyelids, and surrounding tissue (orbicularis oculi muscle);
- Lips (Labium superius oris and labium inferius oris);
- Throat and chest; and
- Any areas that **cannot** sense touch.